

# Compatibility Overview

## Officially Validated Devices

### saremco print CROWNTEC

saremco

CROWNTEC

#### 3D PRINTER

#### CURING UNIT

 ACKURETTA	SOL, DENTIQ, Freeshape 120 G3	3
 ASIGA	MAX UV, PRO 4K	1 2 5
 DentaFab	SEGA	4
 Desktop Health	D4K, Micro Plus XL	1
 DMG	3Delite, 3Demax	1 2 9
 HEYGEARS	UltraCraft Chairside, UltraCraft A2D, UltraCraft A2D 4K	11
 MICROLAY <sup>®</sup>	Versus 385	1
 MIICRAFT	Alpha, Prime	1
 nexa3D	XiP	1 15
 NextDent <small>by 3D SYSTEMS</small>	NextDent 5100, LCD1	7
 phrozen	Sonic 4K 2022, Sonic 4K XL 2022	8
 PRUSA RESEARCH <small>by JOSEF PRUSA</small>	Medical One	12
 rapidshape	D20 II, D30 II, D40 II, D10+, D20+ cartridge, D20+, D30+, D50+	1 2 9
 RAYSHAPE <sup>®</sup>	Shape 1+ Dental	13
 SHINING 3D <sup>®</sup>	AccuFab-L4D/L4K, AccuFab-D1s	6 10
 UNIZ	NBEE	14
 W2P	SF150, SF350, SF650	1

1	<b>OtoFlash by NK-Optik</b> 2 x 2000 Flashes
2	<b>HiLite Power by Kulzer</b> 2 x 180 s
3	<b>CURIE by Ackuretta</b> 2 x 3 min   T:6 P:16 D:10 B:ON <b>CURIE Plus by Ackuretta</b> 2 x 2 min
4	<b>Cure by DentaFab</b> 1 x 10 min
5	<b>Photopol by Dentalfarm</b> 2 x 4 min with vacuum
6	<b>BB Compact by Meccatroncore</b> 1 x 6 min   Power 100%   Temp. 20°C / 68°F   Plate on top rail
7	<b>LC-3DPrint Box by NextDent</b> 1 x 30 min <b>NextDent Cure by NextDent</b> 1 x 10 min
8	<b>Cure by Phrozen</b> 2 x 5 min
9	<b>RS cure by Rapidshape</b> 1 x 6 min
10	<b>Fab Cure by Shining</b> 2 x 10 min <b>Fab Cure 2 by Shining</b> 2 x 5 min (without heater)
11	<b>UltraCraft PCU 3.0 by HeyGears</b> 2 x 3 min 40% power <b>UltraCraft AirCure by HeyGears</b> 2 x 3 min 10% power
12	<b>Medical CW One by Prusa</b> 2 x 10 min
13	<b>ShapeCure by RayShape</b> 2 x 10 min (no heat necessary)
14	<b>U Cure by UNIZ</b> 2 x 2 min at Level 2 Set manually: '+' sign, 'LEVEL', 'CURE Time', 'Save'
15	<b>xCure Desktop by Nexa 3D</b> 2 x 2 min <b>Wash &amp; Cure by Nexa 3D</b> 2 x 10 min

# Good to Know

## saremco print CROWNTEC



### WAYS OF CLEANING

**NEW**



#### RECOMMENDED | EASIEST WAY SAREMCO PRINT CLEANING CONCENTRATE

Mix the concentrate to obtain a cleaning liquid.  
Cleaning in ultrasonic bath (2 x 3 minutes).

More product information about  
**saremco print CLEANING CONCENTRATE**  
can be found on our website.



#### Mixing Protocol



#### Cleaning Protocol



#### OPTIONAL | BEST RESULTS CLEANING BY HAND

Remove excess resin from print object using a **brush and cloths** with a small amount of **IPA (96%)**.

Dry the **surface and inside** of the print object with compressed clean air.

### REACH THE FINAL COLOR

#### RECOMMENDED AUTOCLAVE

**5 minutes Autoclave** (134°C / 276.2°F),  
to accelerate color finalization.



or

#### RECOMMENDED WATER BATH

at least **2 minutes** in **boiling water** (100°C / 212°F),  
to accelerate color finalization.



#### RECOMMENDED POLIMERIZATION LAMP

Expose **each side** of the print object  
for **2 x 20 seconds** at **full power**.



Recommended polymerization devices for light-curing materials such as OtoFlash and HiLite Power reach a wavelength range of 320 - 500nm. Other LED polymerization devices mentioned before, may do not reach the upper wavelength range and do not completely finish the esthetic color finalization process (no influence on the physical properties of the material).

To accelerate the color finalization, it is recommended to place the print object for 5 minutes into an autoclave (134°C / 276.8°F) or in boiling water (100°C / 212°F) for at least 2 minutes. A polymerization lamp such as a Bluephase® G2 from Ivoclar Vivadent can be used optionally (2 x 20 s full power on the print object per side) to cover the higher wavelength ranges (up to 500nm). Please note: This compatibility overview does not replace the instructions for use. Please read the instructions for use carefully.

\*Data Protection Disclaimer

By scanning the QR code, you acknowledge and agree that you will be redirected to YouTube and will be subject to their privacy policy. Please review YouTube's privacy policy before proceeding.