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Instruction for use els composite

1. Product description

Light-curing, highly filled, radio-opaque micro hybrid composite with extremely low shrinkage stress for permanent anterior and posterior restorations. els composite is a polymer-based restorative material type I, class 2, group 1 and 2 according to DIN EN ISO 4049.

2. Composition

Inorganic filler (barium glass and silica 74 % wt, 49 % by volume, particle size between 0,04 - 3,00 µm, median 0,7 µm). BisGMA, BisEMA, catalysts, inhibitors, pigments

3. Intended use

Restorative material intended for the reconstruction or correction of natural dentition.

4. Indication

- Restoration of class I, II, III, IV and V cavities on anterior and posterior teeth
- Extra orally tempered inlays and onlays (direct and indirect), veneers
- Aesthetic corrections of interdental spaces, enamel hypoplasia, discolorations
- Blockings, fractures

5. Contra-indication

Opened pulp, pulpit, known allergy to methacrylates.

6. Patient target group

els composite is suitable for use in all patients without any age or gender restrictions.

7. User

els composite should only be used by a professionally trained dental practitioner.

8. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the pulp with a temporary filling or pulp with suitable underfilling materials (preferably with calcium hydroxide). Note: els composite does not contain TEGDMA or HEMA.

9. Interactions

Avoid underfilling materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

10. Processing stages

10.1. Tooth cleaning

Brush the tooth being treated and the teeth next to it with fluoride-free toothpaste. If necessary, clean interdental spaces using strips and dental floss.

10.2. Shade selection

SAREMCO shades relate to VITA standards.

10.3. Drying

Dry sufficiently, apply rubber dam.

10.4. Cavity preparation

Prepare the cavity as usual. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

10.5. Underfilling

In case of near-to-pulp preparations, it is advisable to protect the pulp using an underfilling (e.g., calcium hydroxide). Do not use zinc oxide eugenol!

10.6. Adhesion

The adhesive has to be conditioned and applied in correspondence to instruction for use of the product. SAREMCO recommends *cmf adhesive system, els duobond* and *els unibond*.

10.7. Restoration with els composite

Application of the syringe: After each usage, turn the syringe screw back by one turn and re-attach the screw cap.

Application of the tip: Remove tip cap. Push the tip into the corresponding recess of the injector and fix and click it into place. (SAREMCO recommends *SAREMCO applicator*). Apply material with constant pressure of the lever.

Gradually adapt *els composite* up to a maximum layer thickness of 2 mm. Light-cure every layer for 40 seconds. When using matrices, remove them and expose filling to light again from lingual and buccal. Note: els composite produces an extremely low shrinkage stress and shrinkage.

Warning: Protect light-curing products from strong sources of light.

The specified exposure time refers to halogen- or LED-light-curing devices with a minimum light intensity of 500 mW/cm² and a wave-length of 400 - 500 nm. They apply to a maximal layer thickness of 2 mm. The necessary exposure time may vary depending on the light source and its instructions for use. In case of doubt, check the lamp's light output and the necessary exposure time before operation in vitro. During polymerization, an inhibition layer forms on the surface; this must not be touched or removed if other composite layers are going to be applied.

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10.8. Finishing, polishing

Prepare the filling with 40 µ and 12 µ diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers. The filling can be prepared and polished immediately after polymerization.

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5. Contre-indication

Close containers after each use with the right lid. Keep out of reach of children. For dental use only. When using tips, make sure that before usage the cap has been removed. Commercially available medical gloves do not provide protection against the sensitization effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice. When polishing or removing composites, it is recommended to always use a water-cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes. Light-curing restorative materials can be processed quickly, as the polymerization can be started by the ambient light of the dental treatment lamp.

6. Précaution

Close containers after each use with the right lid. Rinse thoroughly with water. Consult an eye specialist.

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SAREMCO Dental AG