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### 1. Produktbeschreibung

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#### 8.8. Finishing, Polieren

Flowable, light-curing, radio-opaque microhybrid composite with extremely low shrinkage stress. DIN EN ISO 4049-compliant.

### 1. Product description

Flowable, light-curing, radio-opaque microhybrid composite with extremely low shrinkage stress. DIN EN ISO 4049-compliant.

### 2. Composition

inorganic filler (barium glass and silica 53 % wt, 28 % by volume, particle size between 0.04 - 3.00 µm, median 0.7 µm), BisGMA, BisEMA, catalysts, inhibitors, pigments

### 3. Intended Use

SAREMCO restoratives are intended for the direct reconstruction or correction of functionally compromised natural dentition (e.g., deficient teeth).

### 4. Indication

1. Restorations with minimally invasive preparation technique  
2. Restorations of small cavities and extended fissure sealing

3. Alternate restorations for undercut cavities

4. Restorations of class III - V including wedged shaped defects and cervical caries

5. Repair of fillings, veneers and methacrylate-based temporary restorations

6. First layer of fillings for Class I and II

7. Interlocking of loosened teeth

8. Adhesive attachment of indirect composite and ceramic restorations

9. Adhesive attachment of lingual retainers

### 5. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates

### 6. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is advisable to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

Note: els flow does not contain TEGDMA or HEMA.

### 7. Instructions

Avoid underfilling materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

### 8. Processing stages

#### 8.1. Tooth cleaning

Brush the tooth being treated and the teeth next to it with fluoride-free toothpaste. If necessary, clean interdental spaces using strips and dental floss.

#### 8.2. Shade selection

SAREMCO shades relate to VITA standards.

#### 8.3. Drying

Dry sufficiently, apply rubber dam.

#### 8.4. Cavity preparation

Prepare the cavity as usual. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

#### 8.5. Underfilling

In case of near-to-pulp preparations, it is advisable to protect the pulp using an underfilling (e.g., calcium hydroxide). Do not use zinc oxide eugenol!

#### 8.6. Adhesion

The adhesive has to be conditioned and applied in correspondence to instruction for use of the product. SAREMCO recommends *cmf adhesive system, els duobond and els unibond*.

#### 8.7. Restoration with els flow

**Application of the syringe:** Remove syringe cap and screw on luer-lock cannula. Apply material with constant pressure on the plunger. Do not retract the plunger after application (to prevent a possible slumping of the material).

A retraction of the plunger causes air to be sucked into the material, thereby deteriorating the material properties.  
**Application of the tip:** Remove tip cap. Push the tip into the corresponding recess of the injector and fix and click it into place. (SAREMCO recommends *SAREMCO applicator*). Apply material with constant pressure of the lever.

Gradually adapt els flow up to a maximum layer thickness of 2 mm. Light-cure every layer for 40 seconds. When using metal matrices, remove them and expose filling to light again from lingual and buccal.

**Warning:** After each usage, the application tip must be removed, and the syringe tightly sealed with the original cap. Protect light-curing products from strong sources of light.

The specified exposure times refer to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm<sup>2</sup> and a wavelength of 400 - 500 nm. They apply to a maximal layer thickness of 2 mm. The necessary exposure time may vary depending on the light source and its instructions for use. In case of doubt, check the lamp's light output and the necessary exposure time before operation *in vitro*. During polymerization, an inhibition layer forms on the surface; this must not be touched or removed if other composite layers are going to be applied.

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### 8.8. Finishing, polishing

Prepare the filling with 40 µm and 12 µ diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers. The filling can be prepared and polished immediately after polymerization.

### 9. Storage

Protect light-curing products from strong sources of light and heat! els flow was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F, store in the refrigerator. Do not freeze!

### 10. Batch number and expiry date

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

### 11. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. For dental use only. When using tips, make sure, that before usage the cap has been removed. Commercially available medical gloves do not provide protection against the sensitization effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice. When polishing or removing composites, it is recommended to always use a water-cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

### 12. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

### 13. Hygiene

Use application instruments for one patient only. Do not reuse instruments from patients to avoid contamination. Sterilization of syringes or tips before first use is not required.

### 14. Warranty

The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

### 15. Other notes for Europe

If the user and/or patient become aware of serious incidents connected with application of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides.

### 16. Production / distribution

SAREMCO Dental AG  
Gewerbestrasse 4  
CH-9445 Rebstein / Switzerland  
Tel: +41 (0) 71 775 80 90  
Fax: +41 (0) 71 775 80 99  
info@saremco.ch  
www.saremco.ch

Edited 02/2022 | D600166

Class IIa medical devices  
Class III medical devices (Canada)

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### 2. Composition

inorganic filler (barium glass and silica 53 % wt, 28 % by volume, particle size between 0.04 - 3.00 µm, median 0.7 µm), BisGMA, BisEMA, catalysts, inhibiteurs, pigments

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