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### 1. Produktbeschreibung

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## Instruction for use els duobond

### 1. Product description

Dual-curing two-component self-etching bond to create a permanent marginal-gap-free adhesion between the tooth structure (dentine, enamel) and the light-, self- or dual-curing filling/fixing material.

### 2. Composition

els duobond base: BisEMA, methacrylated phosphoric salt, initiators  
els duobond catalyst: acetone, water, initiators

### 3. Indication

- Direct self- or dual-curing composite restorations and core build-ups
- Direct restorations with light-curing composite-based materials
- Fixing root pins with dual- or self-curing composite cements
- Indirect restorations: at the use of dual-curing or self-curing composite cements to fix inlays, onlays, crowns and bridges.

### 4. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates.

### 5. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide). Note: els duobond does not contain TEGDMA or HEMA.

### 6. Interactions

Avoid underfilling materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol, belong to this category.

### 7. Processing stages

#### 7.1. Drying

Due to the following application of the filling material, a drainage of the working field is required. Avoid any contamination of the cavity with blood or saliva. A rubber dam is recommended.

#### 7.2. Cavity preparation

Prepare the cavity as usual. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity. Avoid overdrilling.

#### 7.3. Mixing

els duobond is dual-curing and therefore must be processed immediately after mixing. Thoroughly mix 1 drop of the base and 1 drop of the catalyst on a mixing tray with a brush for around 5 seconds (forms a self-etching adhesive). After use, rub the adhesive seal with alcohol to avoid contamination.

#### 7.4. Application

Apply els duobond to enamel/dentine in a not too thin layer with a brush and massage into the tooth substance for 20 seconds. Note: els duobond is light-curing, so avoid too intense ambient light. Dim the surgical light during the application.

#### 7.5. Drying

Dry the adhesive layer with an air blower for at least 5 seconds.

#### 7.6. Polymerizing

##### 7.6.1. Direct restorations

When using light-, self- or dual-curing filling material polymerize the bonding layer for 10 seconds with blue light (halogen or LED light).

##### 7.6.2. Indirect restorations

**Light-curing luting composites:** When using light-curing luting composites polymerize the bonding layer for 10 seconds with light (halogen or LED light).  
**Self- or dual-curing luting composites:** When using dual- or self-curing luting materials, els duobond creates good adhesion even without light-curing. Optionally, light-curing can be additionally performed for 10 seconds with blue light (halogen or LED lamp).

### 7.7. Fixing root pins

#### 7.7.1. Bonding

To fix root pins, the mixed els duobond (base and catalyst) is massaged into the root canal and onto the occlusal surface of the stump with an applicator for 20 seconds. Disperse solvent with oil-free air for 5 seconds. Remove any excess material by dabbing with an absorbent paper strip. Do not light-cure the bonding layer!

#### 7.7.2. Inserting

Insert the root pin prepared according to the manufacturer's instructions into the root canal using a suitable fixing material in accordance with the instructions for use and to achieve slight excess cement. SAREMCO recommends els cem.

#### 7.7.3. Light-curing

Light-cure for at least 40 seconds to fix the root pin and to harden the occlusal applied els duobond.

### 7.8. Notes for fixing Maryland bridges

When fixing Maryland bridges with els duobond, the enamel must be roughened or etched with phosphoric acid gel (e.g. SAREMCO cmf etch) prior to use of the bond. Self-etching bonds are less effective on non-prepared enamel. Excessive amount of bond material onto unpolished and unbevelled enamel can cause edge staining. For optimal bonding results with els duobond, the enamel should be prepared in advance in class III, IV, V and for diastema closures.

### 8. Storage

Close els duobond bottles tightly after use. Store in an upright position to ensure the backflow of the liquid. Do not expose light-curing products to direct sunlight or operating light. els duobond was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store in a refrigerator at 4°C - 8°C / 39°F - 46°F. Do not freeze! Bring to room temperature before use. Constant temperatures above 25°C / 77°F can reduce the shelf-life of the product.

### 9. Batch number and expiry date

The batch number should be specified to identify products in the case of enquiries. Products should no longer be used once the expiry date has elapsed.

### 10. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes into contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice.

### 11. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

### 12. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination.

### 13. Warranty

Our liability is restricted to the quality of our products. In case of a product being of defective quality, only its value is replaced. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental.

### 14. Other notes for Europe

If the user and/or patient become aware of serious incidents connected with application of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides.

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