Gebruichsanweisung
saremco print - CROWNTEC

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1. Product description
CROWNTEC is a light-curing, flowable polymer based on methacrylic acid ester for production of 3D-printed permanent crowns, inlays, onlays and veneers, temporary crowns and bridges and artificial teeth.

2. Composition
SAREMCO dental glass (silanized), pyrogenic silica, catalysts, inhibitors.

3. Indication
With the aid of the 3D-printing machine from ASIGA and Rapid Shape:
1. Production of permanent crowns, inlays, onlays and veneers
2. Production of temporary crowns and bridges, inlays, onlays and veneers
3. Production of artificial teeth for subsequent insertion into a denture base


4. Contra-indication
Do not use the product in case of a known allergy to one or more ingredients. In case of doubt, clarify and exclude a possible allergy with the help of a specific allergy test before using CROWNTEC.

5. Side effects
In individual cases, contact allergies to products with similar composition have been described.

6. Interactions
None known.

7. Processing stages
The following instructions apply to the model modelled on the computer:

Minimum wall thickness: The following illustration shows the specified minimum wall thicknesses for the respective indication: the wall thickness must not be undercut even after manual grinding.

The following applies to temporary bridges: connector area at least 16 mm².

7.1. Generating printing file
Generate the printing file of the desired restoration by using appropriate software (Composer) and deliver it suitable to the printer.

Important Note: commercial available artificial teeth may be subject of copyright law. When using a st-file of such teeth, copyright law must be considered.

7.2. Printing
Work as clean as possible, as dirty reservoirs or machines can cause deformation/discolouration and therefore failure of the printed objects.

7.3. Cleaning
After the printing process is completed, remove the building platform from the machine.

7.4. Finishing the printed jobs
To achieve the desired material properties and biocompatibility, post-curing of the completely dried and cleaned printed objects is necessary.

Warning: Protect light-curing products from strong light sources!

7.5. Fastening the definitive crowns, inlays, onlays and veneers
In case of definitive single crowns, the inside of the crowns should be roughened with a sandblast (Al2O3, 110 µm). Then, as usual, fix it definitively with a composite cement material.

7.5.2. Fastening the temporary crowns and bridges, inlays, onlays and veneers
Fasten the finished transitional prosthesis with commercially available provisional cements.

7.5.3. Connecting of artificial teeth and prosthesis
a) Inserting the printed artificial teeth in a printed, prefabricated denture base. Roughen the base surface of the printed artificial teeth for example by sandblasting (Al2O3, 110 µm), apply a primer and a fixing material, insert in the prosthesis according to the natural shape and polymerise.

7.6. Finishing, polishing
Prepare the restoration with 40 µm and 12 µm diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers.

8. Storage
Protect this product from strong light and heat sources! The recommended storage temperature is between 4°C and 28°C / 39°F and 82°F.

9. Batch number and expiry date
The batch number is used to identify the product in case of queries.

10. Precautionary measures
For dental use only. Keep out of reach of children. Use of nitrile gloves while working with CROWNTEC is recommended.

11. Emergency measures
In case of direct contact of the uncured material with the oral mucosa, rinse with water.

12. Hygiene
Restorations made of CROWNTEC should not be decontaminated with chemical products.

13. Warranty
Our liability is restricted to the quality of our products. In case of a product being of defective quality, only its value is replaced.

14. Other notes for Europe
If the user and/or patient become aware of serious incidents connected with application of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides.

15. Production / distribution
SAREMCO Dental AG
Gewerbestrasse 4
CH-9445 Rebstein / Switzerland

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