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## Instruction for use saremcoprint - CROWNTEC

**1. Product description**  
CROWNTEC is a light-curing, flowable polymer based on methacrylic acid ester for production of 3D-printed permanent crowns, inlays, onlays and veneers, temporary crowns and bridges and artificial teeth.

**2. Composition**  
BisEMA, dental glass (silanized), pyrogenic silica, catalysts, inhibitors.

**3. Intended use**  
saremcoprint products provide light-curing 3D printable resin-based materials for the correction or reconstruction of functionally compromised natural dentition (e.g., missing teeth or deficient teeth) by manufacturing of customized 3D-printed dental prostheses.

**4. Indication**  
With the aid of the 3D-printing machine from ASIGA and Rapid Shape:

- Production of permanent crowns, inlays, onlays and veneers
- Production of temporary crowns and bridges, inlays, onlays and veneers
- Production of artificial teeth for subsequent insertion into a denture base

**5. Contra-indication**  
Do not use the product in case of a known allergy to one or more ingredients. In case of doubt, clarify and exclude a possible allergy with the help of a specific allergy test before using CROWNTEC.  
CROWNTEC must not be used for any other purposes than those specified in the "Indication" section. Any deviation from this instruction for use may have negative effects on the chemical and physical quality of the restorations produced from CROWNTEC.

**6. Patient target group**  
CROWNTEC is suitable for use in all patients without any age or gender restrictions.


**7. User**  
CROWNTEC should only be used by a professionally trained dental practitioner.

**8. Side effects**  
In individual cases, contact allergies to products with similar composition have been described.

**9. Interactions**  
None known.

**10. Processing stages**  
Please also refer to the Fabrication Manual for ASIGA and Rapid Shape at [www.saremc.ch](http://www.saremc.ch).

The following instructions apply to the model modelled on the computer:  
Minimum wall thickness: The following illustration shows the specified minimum wall thickness for the respective indication: the wall thickness must not be undercut even after manual grinding.  
The following applies to temporary bridges:  
connector area at least 16 mm<sup>2</sup>.  
The connector area should be at least as large as possible.  
For physical stability, the height of the connector is more important than the width. Doubling the width results in only doubling the strength, while doubling the height results in eight times the strength. Oval connector faces are recommended.

Design temporary (long-term) bridges in the non-visible molar region in the form of a floating bridge (posterior bridge). A floating bridge does not sit firmly on the jaw, but forms a surface that can be rinsed underneath and can therefore be optimally cleaned

**10.1. Generating printing file**  
Generate the printing file of the desired restoration by using appropriate software (Composer) and deliver it suitable to the printer. Please observe the corresponding instruction for use of software and printer.

**Important note:** artificial teeth that are commercially available may be subject to copyright law. When using a print file of those teeth, copyright laws must be considered.

**10.2. Printing**  
Work as clean as possible, as dirty reservoirs or machines can cause deformation/discolouration and therefore failure of the printed objects.  
Briefly shake the liquid material and pour it into the reservoir of the 3D-printing machine. Start the printing process by following the instruction for use of the printer. Download the parameter set intended for CROWNTEC from the database of the printer manufacturer. For ASIGA printers, the following applies: A working temperature of 35°C/95°F must be maintained. For all other printers, it is recommended that both the printer and the resin be brought to operating temperature. A cold start should be avoided.

For saremcoprint compatibility overview of officially validated devices, scan QR-Code using camera app or use scan apps such as Lens (Google) or QR-Code-Scanner (QR SCAN Team). Download the compatibility overview (use Mozilla Firefox or Google Chrome) from the homepage [www.saremc.ch/en/download/instructions-for-use/](http://www.saremc.ch/en/download/instructions-for-use/). Look under compatibility overview.

**10.3. Cleaning**  
After the printing process is completed, remove the building platform from the machine. During removing the restoration and the following cleaning steps, wearing glove (nitrile gloves) and protective goggles are advised.  
Place the platform on a piece of paper or cloth with the built jobs facing upwards. Remove the printed jobs from the platform by using a suitable instrument (putty knife). The support structures are then cut off. Either a cut-off wheel or a side cutter can be used to do this. To remove excess material, clean the printed job with an alcohol-soaked (96%) cloth and possibly a brush soaked in an alcohol solution until all resin remains are completely removed. Then dry the printed jobs thoroughly with an air syringe.

**Warning:** Protect light-curing products from strong light sources!

**10.4. Finishing the printed jobs**  
Optional step 1: Carefully blast the surface of the print objects with a blast polishing material. Wearing nitrile gloves, protective goggles, and a dust mask is recommended for the finishing step.

Optional step 2: The objects can be individualized by means of stains, e.g., els paintart. Observe the instructions for use.

To achieve the desired material properties and biocompatibility, post-curing of the completely dried and cleaned printed objects is necessary. For final polymerization place the printed jobs in a UV-light box.  
Note: time of curing depends greatly on type of lamps / lightbox used. The final properties and the final colour depend on the post-curing process. Post-curing is an UV-light treatment to ensure that saremcoprint materials obtained full polymer conversion, the residual monomer is reduced to a minimum and the highest mechanical properties are achieved.

This procedure is a necessary step to attain a biocompatible end-product. It is suggested to use the polymerization units "Signum HiLite Power" from Heraeus Kulzer (2 x 180s) or the UV-Flash device "Otoflash G171" from NK-Optik (4000 flashes). In general, all lightboxes for light-curing veneering materials can be used that cover a wavelength range of 320 - 500 nm.  
Lightboxes with integrated flash light allow shorter exposure time compared to conventional lamps. Always follow the respective instruction of use of the polymerization units. The times indicated refer to light curing units that are regularly maintained and tested for light intensity.

**10.5. Fastening**  
**10.5.1. Fastening the definitive crowns, inlays, onlays and veneers**  
In case of definitive single crowns, the inside of the crowns should be roughened with a sandblast (Al<sub>2</sub>O<sub>3</sub>, 110 µm). Then, as usual, fix it definitively with a composite cement material. Zinc-phosphate cements as well as glass-ionomer-cements are only of limited suitability due to their opacity. The fastening composites Panavia 5 [Kuraray] and Variolink [voclar] are recommended.

**10.5.2. Fastening the temporary crowns and bridges, inlays, onlays and veneers**  
Fasten the finished transitional prosthesis with commercially available provisional cements.

## EN

**10.5.3. Connecting of artificial teeth and prosthesis**  
a) Inserting the printed artificial teeth in a printed, prefabricated denture base. Roughen the base surface of the printed artificial teeth for example by sandblasting (Al<sub>2</sub>O<sub>3</sub>, 110 µm), apply a primer and a fixing material, insert in the prosthesis according to the natural shape and polymerise.  
Alternatively, CROWNTEC can also be used directly as fixing material. Therefore put a small amount of material with a brush on the roughened teeth-surface of the artificial tooth, put it into the prosthesis, eliminate any excess material and light cure it from all sides for at least 20 seconds. The polymerization light should have at least a light output of 600 mW/cm<sup>2</sup>.  
b) Using a classical finishing procedure like the pouring method with cold cure resin after roughening the teeth.

**10.6. Finishing, polishing**  
Prepare the restoration with 40 µ and 12 µ diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers.

**11. Storage**  
Protect this product from strong light and heat sources! The recommended storage temperature is between 4°C and 28°C / 39°F and 82°F. Close the package after each use.

**12. Batch number and expiry date**  
The batch number is used to identify the product in case of queries. Do not use this product after the expiration date.

**13. Precautionary measures**  
For dental use only. Keep out of reach of children. The use of nitrile gloves while working with CROWNTEC is recommended until post-curing. Commercially available medical gloves do not provide effective protection against the sensitising effect of methacrylates.  
If the product comes into contact with the glove, remove the glove and dispose of it, wash your hands immediately with soap and water and put on a new glove. In case of an allergic reaction, consult a doctor. When polishing or removing composites, it is recommended to always use a water cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

**14. Emergency measures**  
In case of direct contact of the uncured material with the oral mucosa, rinse with water. In case of contact with eyes, rinse thoroughly with water and consult an eye specialist.

**15. Hygiene**  
Restorations made of CROWNTEC should not be cleaned with chemical products. Cleaning with water is sufficient. The finished restorations can - if necessary - be disinfected with an ethanol solution.

**16. Warranty**  
Our liability is restricted to the quality of our products. In case of a product being of defective quality, only its value is replaced. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental.

The Summary of Safety and Clinical Performance of SAREMCO Products can be found in the European database on medical devices (EUDAMED-<https://ec.europa.eu/tools/eudamed>).

**17. Production / distribution**  
SAREMCO Dental AG  
Gewerbestrasse 4  
CH-9445 Rebstein / Switzerland  
Tel: +41 (0) 71 775 80 90  
Fax: +41 (0) 71 775 80 99  
info@saremc.ch  
www.saremc.ch

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