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## Instruction for use zirconia primer

### 1. Product description

Adhesive for zirconium dioxide frameworks. *zirconia primer* ensures a strong bonding with acrylates and zirconia ceramics and is designed for use in the dental laboratory for indirectly manufactured dentures.

### 2. Composition

Ethanol, MDP, Methacryloxypropyltrimethoxysilan

### 3. Intended Use

Saremco cements and auxiliary bonding promoters are intended for the indirect reconstruction or correction of functionally compromised natural dentition (e.g. deficient teeth).

### 4. Indication

Priming of frameworks of high-performance ceramics such as zirconium dioxide.

### 5. Contra-indication

Known allergy to methacrylates. *zirconia primer* is not suitable for conditioning composite frameworks.

### 6. Patient target group

*zirconia primer* is suitable for use in all patients without any age or gender restrictions.

### 7. User

*zirconia primer* should only be used by a professionally trained dental practitioner.

### 8. Side effects

In individual cases, contact allergies have been described to products of a similar composition.

### 9. Interactions

Avoid contact with materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

### 10. Processing stages

#### 10.1. Conditioning

Blast the zirconium dioxide frameworks/surfaces with 110 µ aluminum oxide at 2 - 3 bar (roughening, increasing the surface) and clean with oil-free compressed air. Then use a brush to apply the *zirconia primer* to the surface to be conditioned. Dry with oil-free compressed air.

**Caution:** If the air flow is too strong, the primer will be scattered, which will lead to a worse adhesion.

#### 10.2. Further processing

Apply the veneering material/fixing cement to the surface to be treated and cure according to the manufacturer's instructions. Saremco recommends elis cem. Using dual-curing materials, additional light-curing is recommended (also through the respective restoration).

#### 11. Storage

Protect light-curing products from strong sources of light and heat! *zirconia primer* was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Constant temperatures above 28°C / 82°F can reduce the shelf-life of the product.

#### 12. Batch number and expiry date

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

#### 13. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. For dental use only. Commercially available medical gloves do not provide protection against the sensitization effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice.

**Hazard warnings:** May cause an allergic skin reaction. Causes serious eye damage. Harmful to aquatic organisms, with long-term effect.

#### 14. Emergency measures

In case of contact with the eyes, rinse cautiously with water for a few minutes. Remove any contact lenses, if possible. Continue rinsing. Consult an eye specialist.

#### 15. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination. The usual hygiene rules must be observed. Dispose of the content/container in accordance with local, regional, national, and international regulations.

#### 16. Warranty

The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user

to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of Saremco Dental.

### 17. Other notes for Europe

If the user and/or patient becomes aware of serious incidents in connection with the use of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides.

The Summary of Safety and Clinical Performance of Saremco Products can be found in the European database on medical devices (EUDAMED-<https://ec.europa.eu/tools/eudamed>).

### 18. Production / distribution

SAREMCO Dental AG  
Gewerbestrasse 4  
CH-9445 Rebstein / Switzerland  
Tel: +41 (0) 71 775 80 90  
Fax: +41 (0) 71 775 80 99  
[info@saremco.ch](mailto:info@saremco.ch)  
[www.saremco.ch](http://www.saremco.ch)

Edited 10-2022 | D600234

Class IIa medical devices  
Class III medical devices (Canada)

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