



## apt<sup>®</sup> advanced polymer technology unibond

### Technical Data

Resin	HEMA- and TEGDMA-free, lightcuring 1-Component Self-Etch-Adhesive.
Composition base	ethanol, methacrylated phosphoric salt, aromatic urethanmethacrylate, BisEMA, water, fillers, initiators
Consistency (mixed)	low viscous
Shear bond strength dentine	> 10 MPa (bovine dentine)
Shear bond strength enamel	> 10 MPa (bovine enamel)

Version: 12-2019