

csp ceramic silane

Technical Data

Product description	Two components primer for the silanisation of silicatic or aluminiumoxidic ceramics
Composition	
silane base:	gamma-methacryloxypropyltrimethoxysilane
silane activator:	water, acetic acid
Mixing ratio base : activator	2 : 1
Consistency	low viscous
Shearbond strength on polished untreated feldspathic ceramic (24h) ¹⁾	> 8 MPa

1) M. Cattani-Lorente, A. Wiskot, ADHESIVE EFFICACY OF TWO SILANES, Report 2010, University of Geneva, SECTION DE MEDECINE DENTAIRE, DEPARTEMENT DE PROTHESE DENTAIRE

Version: 12-20219