

james-2 bond

Technical Data

Definition	Solvent-free, light curing single component adhesive (primer and bond) for enamel and dentine
Composition	Hydroxyethylmethacrylate Urethanemethacrylate Polyalkenoate methacrylized Hydroxypropylmethacrylate Glycerinedimetharylate
Appearance	Clear yellowish liquid
Consistency	Viscous liquid
Refractive index (25°C)	1.458-1.468
Shear Bond Strength	16.9 MPa (Degrange, 2003)

Version: 12-2019