

**saremco print DENTURETEC (provisional)**

**Technical Data**

|                                     |                                  |
|-------------------------------------|----------------------------------|
| Consistency                         | middle viscos homogeneous liquid |
| Colour                              | Gum-coloured                     |
| Depth of cure (DIN EN ISO 4049)     | ≥ 1,5 mm                         |
| Flexural strength (DIN EN ISO 4049) | > 100 MPa                        |
| E-Modulus (DIN EN ISO 4049)         | > 1000 MPa                       |
| Barcol-Hardness                     | > 20                             |
| Shelf-life                          | 2 years                          |

Version: 05-2019